

# MEDIPO - ZT s.r.o.

ASSOCIATION OF MANUFACTURERS AND SUPPLIERS OF MEDICAL DEVICES

office I.: Kohoutova 24, 613 00 Brno  
tel./fax: +420 548 525 348  
e-mail: medipo-zt@volny.cz

office II.: Slaměnikova 23, 614 00 Brno  
tel.: +420 545 210 119  
web: www.medipo.cz

The company MEDIPO – ZT develops, produces and sales a manufactures of medical devices. It means the technical equipment for infusion, lythotripts and endoscopes.



The infusion stand with two brakes.

## INFUSION PUMP IT 1013


The IT 1013 infusion pump is intended for intravenous and intra-arterial infusion of medicaments, parenteral or enteral nutrition, or blood transfusion. A correct run of the infusion is controlled by dropping sensor, air bubble detector and pressure sensor in the hose of the infusion set. The IT 1013 may be attached to a mobile stand, vertical or horizontal bar, or a stand on a rigid support plate (shelf, table, etc.)





## CHILDREN'S VAGINOSCOPE MT 1000

The MT 1000 children's vaginoscope is designed for an easy gynecological examination of children. With the vaginoscope it is possible to see changes in vaginal mucosa, cervix uteri, or location of foreign object. The vaginoscope consists of a set of tubes of various diameters and lengths, equipped with appropriate obturators and optics. Illumination is provided by halogen light source with the light-guide cable.



## REKTOSKOP MT 1100

The MT 1100 rectoscope is calibrated rigid endoscope, designated for examination, extraction of tissue samples and therapeutic intervention in the rectum area. The light is fed into the rectoscope from the source of cool light by the means of flexible light-guide cable.



## MEDILIT M

The lithotripter MEDILIT is a powerful equipment, intended for non-invasive disintegration of kidney and gall-bladder stones and stones in the urinary and biliary tract using focused shock-wave. Whole process of the therapy is computer-controlled. There are two localization systems, both X-ray and ultrasound.